

# MicroSnap™

A new light on Pathogen Detection



MicroSnap™ is a rapid bioluminogenic test method for the detection and enumeration of Coliform and E.coli bacteria. The new MicroSnap™ platform replaces traditional microbiology testing methods with a rapid, scientific test that provides results in less than eight hours. Same-shift results allow food and beverage processors to screen materials faster, monitor the plant environment in real time and release finished products faster. The MicroSnap™ consists of an Enrichment Device containing a growth medium and a Detection Device containing a bioluminogenic (light-producing) substrate. When MicroSnap™ detects the specific microorganism, light is emitted and measured with the EnSure™ luminometer.



## Features & Benefits:

- Fastest microorganism test for Coliform and E. coli available
- Allows for faster quality assurance response and product release.
- Quantitative (enumeration) results in **6 hours**
- Qualitative (presence/absence <1CFU/100mL) results in **8 hours**
- Simple pass/fail result at desired specification
- Low level detection (<10 organisms)
- No special sample preparation required
- Self-contained devices provide ease of use
- Unique liquid-stable reagent provides high sensitivity and repeatability
- Uses proven conventional diagnostic properties
- Independent of sample effects
- Simple 2 step procedure
- Snap-Valve™ technology & Built-In-Pipette provide ease of use
- Equivalent results to other cultural methods
- Repeatable testing on same enrichment sample

## Applications Include:

- Surface swabbing
- Raw material and finished product testing
- Solids, liquids and filterable products

## Detection Principle:

**Compound x – Co-Factors + Enzyme C + ATP + Luciferase → Light**




Phone: 416-261-4865  
Email: [quotes@scigiene.com](mailto:quotes@scigiene.com)  
[www.scigiene.com](http://www.scigiene.com)


# Step: 1 Enrichment

MicroSnap™


## Step 1: Environmental Surface Swabs, Liquids and Solid Samples




**1.1 Surface:** Swab a 10x10cm area or larger depending on protocol with the [MicroSnap Enrichment Swab](#) (ATP-ES100)



**1.2 Liquids\*:** 1mL beverage or water sample added directly to [MicroSnap Enrichment Swab](#).




**1.3 Solid Samples:** 1mL 10% w/v suspension of solid samples added directly to [MicroSnap Enrichment Swab](#).



**2.** Reinsert Snap-Valve bulb into swab tube.




**3. Activate the device.** Bend bulb forward and backward to break the Snap-Valve.



**4.** Lift the bulb up (about 1-2") and squeeze the bulb to release the liquid into tube. Release pressure from the bulb (the bulb is like a dropper bulb) and replace bulb in the tube. Most liquid should be in the bottom of the tube.



**5.** Shake the tube gently to mix sample in the liquid.



**6.** Incubate at 37° ±0.5°C in [Dry Block Incubator](#) for 6 hours for a quantitative measurement or 8 hours for a qualitative measurement. This is the Enriched sample. **Proceed to step 2.**


# Step: 2 Detection

MicroSnap™


## Step 2: Detection



**1.** Shake [MicroSnap Coliform test](#) (ATP-CC100), detection tube by tapping on the palm of your hand 5 times to bring liquid to the bottom of the tube.




**2:** Aseptically transfer .1mL (3 drops or to full line) of Enriched sample from [MicroSnap Enrichment Swab](#) to [MicroSnap Coliform Test](#) (ATP-CC100)




**3.** Activate [MicroSnap Coliform Test](#) by bending the bulb forward and backward, breaking the Snap-Valve.



**4.** Shake the tube gently to mix sample in the liquid.



**5.** Incubate [MicroSnap Coliform Test](#) for 10±0.2 minutes at 37° ±0.5°C in [Dry Block Incubator](#)



**6.** Insert [MicroSnap Coliform Test](#) in the [EnSURE](#) luminometer (ATP-206) and initiate the measurement. Record the results as RLUs and refer to Table 1 to interpret results.



**7.** If a positive result is obtained for Coliform, the presence of E.coli can be verified using the [MicroSnap E.coli Test](#) (ATP-EC100). Repeat the Step 2 procedure using another aliquot sample from the same Enriched sample and the [MicroSnap E.coli test](#).

**TABLE 1: RELATIONSHIP BETWEEN ESTIMATED CFU AND MICROSNAP COLIFORM/E. COLI RLU**


Estimated CFU	Equivalent RLU for 10 minutes assay on	
	SystemSURE Plus	EnSURE
<10	2	2
<20	3	4
<50	6	7
<100	8	12
<200	12	20
<500	25	35
<1000	50	60
<5000	85	180
<10,000	150	300



## Results & Interpretation

MicroSnap™

This graph represents the limit of detection of MicroSnap. A range of ready-to-eat foods were inoculated with 500 cfu/ml of *Escherichia coli* and tested at given intervals. Results show that significant levels of *E.coli* were detected after only 5 hours.



Method and Time to Results	<a href="#">MicroSnap</a> 10 cfu.g or ml <b>7 hours + 10 minutes</b>	Traditional Method (VRBG) ISO 16140 <b>24 hours</b>
	% positive	% positive
50 food types from 5 food groups, (Meat, RTE, Salads, Milk and Dried Food)	99	95


MicroSnap is more accurate than traditional methods and yields results in significantly less time.

	Coliforms <a href="#">MicroSnap, ATP-CC100</a>	E.coli <a href="#">MicroSnap, ATP-EC100</a>
<b>Sensitivity/NPV (%)</b>	94	89
<b>Specificity/PPV (%)</b>	99	100

Third party validation shows sensitivity and specificity of the MicroSnap devices is excellent and better than traditional and convenience microbiological methods. **Sensitivity** can be defined as the ability of a method to detect a target organism compared to the ability of the reference method. **Specificity** can be defined as the ability of the method to detect only the target organism, and not suffer interference from non-target organisms, compare to the reference method.

	Testing Method for <a href="#">E.coli</a> detection		
	Petrifilm EC	DryCult Coli	<a href="#">MicroSnap</a>
<b>Sensitivity (%)</b>	40	33	89
<b>Specificity (%)</b>	91	79	100





At very low level contamination Micro-Snap clearly detects the organisms in food matrices.

## Summary

The MicroSnap™ has been proven to:

- ✓ Allow for faster quality assurance and product release
- ✓ Detect Low levels of specific organisms (<10 organisms)
- ✓ Provide quantitative results in **6 hours** and qualitative results in **8 hours**
- ✓ Allow repeat testing on one enrichment sample
- ✓ Be the fastest microorganism test for Coliform and E.coli today!

Best of all, the MicroSnap™ Coliform & E.coli tests are now:



# Product Order Information

Part Number	Description
<a href="#">ATP-ES100</a>	MicroSnap™ STEP 1: Enrichment Swab- E.coli & Coliform (100 count)
<a href="#">ATP-CC100</a>	MicroSnap™ STEP 2: Detection Device for Coliform (100 count)
<a href="#">ATP-EC100</a>	MicroSnap™ STEP 2: Detection Device for E.coli (100 count)
<a href="#">ATP-206</a>	EnSURE™ Luminometer
<a href="#">ATP-625</a>	Small Format Digital Dry Block Incubator with a 12-Swab Well Heating Block
<a href="#">ATP-MS190</a>	<b>*new*</b> Positive Control Kit for MicroSnap Coliform and/or MicroSnap E. coli, 10 vials
<a href="#">PT1A</a>	Timer & Clock -4 event

**To place an order call or email us today:**



Phone: 416-261-4865  
 Email: [quotes@scigiene.com](mailto:quotes@scigiene.com)  
[www.scigiene.com](http://www.scigiene.com)